

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Micro Ortho Plus

Print date: 19.04.2017

Product code: 10-500\_EN

Page 1 of 6

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

**1.1. Product identifier** International Nomenclature Code **39 06 1000**

Micro Ortho Plus

Product group: Powder component

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Creation of dental products.

### 1.3. Details of the supplier of the safety data sheet

Company name: Micropolymer S.r.l.s.  
Via Torrette Battifoglia Zona  
Street: Industriale Fraz.S.Andrea Delle Fratte  
Place: 06132 Perugia (Italy)  
Telephone: +39 3274572561  
e-mail: info@micropolymerdental.com  
Contact person: Ing.Libero Leone  
e-mail: info@micropolymerdental.com  
Internet: www.micropolymerdental.com

### 1.4. Emergency telephone number:

+390266101029

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### Regulation (EC) No. 1272/2008

This mixture is not classified as hazardous according to Regulation (EC) No. 1272/2008.

#### 2.2. Label elements

#### 2.3. Other hazards

Contains: Methyl methacrylate Organic peroxides. May produce an allergic reaction.

### SECTION 3: Composition/information on ingredients

#### 3.2. Mixtures

##### Chemical characterization

Mixture on base of acrylic powder and catalyst.

##### Hazardous components

CAS No	Chemical name	Index No	REACH No	Quantity
	EC No			
	Classification according to Regulation (EC) No. 1272/2008 [CLP]			
	Acrylic polymer			95 - < 100 %

Full text of H and EUH statements: see section 16.

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### General information

Remove contaminated, saturated clothing immediately. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Micro Ortho Plus

Print date: 19.04.2017

Product code: 10-500\_EN

Page 2 of 6

#### After inhalation

Provide fresh air.

#### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

#### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

#### After ingestion

Rinse mouth immediately and drink plenty of water.

#### **4.2. Most important symptoms and effects, both acute and delayed**

No information available.

#### **4.3. Indication of any immediate medical attention and special treatment needed**

Treat symptomatically.

### SECTION 5: Firefighting measures

#### **5.1. Extinguishing media**

##### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

##### **Unsuitable extinguishing media**

Full water jet

#### **5.2. Special hazards arising from the substance or mixture**

Non-flammable.

#### **5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

#### **Additional information**

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

### SECTION 6: Accidental release measures

#### **6.1. Personal precautions, protective equipment and emergency procedures**

Avoid generation of dust. Do not breathe dust.

#### **6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

#### **6.3. Methods and material for containment and cleaning up**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

#### **6.4. Reference to other sections**

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

### SECTION 7: Handling and storage

#### **7.1. Precautions for safe handling**

##### **Advice on safe handling**

No special measures are necessary.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

##### **Further information on handling**

No information available.

#### **7.2. Conditions for safe storage, including any incompatibilities**

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

Print date: 19.04.2017

**Micro Ortho Plus**  
Product code: 10-500\_EN

Page 3 of 6

## Requirements for storage rooms and vessels

Keep container tightly closed.

## Advice on storage compatibility

No special measures are necessary.

## Further information on storage conditions

No information available.

## 7.3. Specific end use(s)

Dental prothesis.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

### 8.2. Exposure controls

#### Protective and hygiene measures

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

#### Eye/face protection

Wear eye/face protection.

#### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

#### Skin protection

Wear suitable protective clothing.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state:	Powder
Colour:	See packaging
Odour:	characteristic

#### Test method

pH-Value:	not determined
-----------	----------------

#### Changes in the physical state

Melting point:	ca. 110 °C
----------------	------------

Initial boiling point and boiling range:	not determined
--	----------------

Flash point:	> 250 °C
--------------	----------

#### Flammability

Solid:	not determined
--------	----------------

Gas:	not applicable
------	----------------

Lower explosion limits:	not determined
-------------------------	----------------

Upper explosion limits:	not determined
-------------------------	----------------

Ignition temperature:	> 400 °C
-----------------------	----------

#### Auto-ignition temperature

Solid:	not determined
--------	----------------

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Micro Ortho Plus

Print date: 19.04.2017

Product code: 10-500\_EN

Page 4 of 6

Gas:

not applicable

Decomposition temperature:

not determined

#### **Oxidizing properties**

Not oxidizing.

Vapour pressure:

not determined

Density (at 20 °C):

1,16 g/cm<sup>3</sup>

Bulk density (at 20 °C):

700 - 750 kg/m<sup>3</sup>

Water solubility:

insoluble

#### **Solubility in other solvents**

Soluble in hydrocarbons.

Partition coefficient:

not determined

Vapour density:

not determined

Evaporation rate:

not determined

#### **9.2. Other information**

Solid content:

100 %

### SECTION 10: Stability and reactivity

#### **10.1. Reactivity**

No hazardous reaction when handled and stored according to provisions.

#### **10.2. Chemical stability**

The product is stable under storage at normal ambient temperatures.

#### **10.3. Possibility of hazardous reactions**

No known hazardous reactions.

#### **10.4. Conditions to avoid**

none

#### **10.5. Incompatible materials**

No information available.

#### **10.6. Hazardous decomposition products**

Exothermal decomposition with formation of: Methyl methacrylate

### SECTION 11: Toxicological information

#### **11.1. Information on toxicological effects**

##### **Toxicocinetics, metabolism and distribution**

The product has not been tested. .

##### **Acute toxicity**

Based on available data, the classification criteria are not met.

##### **Irritation and corrosivity**

Based on available data, the classification criteria are not met.

##### **Sensitising effects**

Based on available data, the classification criteria are not met.

##### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

##### **STOT-single exposure**

Based on available data, the classification criteria are not met.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Micro Ortho Plus

Print date: 19.04.2017

Product code: 10-500\_EN

Page 5 of 6

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

This mixture is classified as not hazardous according to 1999/45/EC.

### SECTION 12: Ecological information

#### 12.1. Toxicity

The product is not: Ecotoxic.

#### 12.2. Persistence and degradability

The product has not been tested.

#### 12.3. Bioaccumulative potential

The product has not been tested.

#### 12.4. Mobility in soil

The product has not been tested.

#### 12.5. Results of PBT and vPvB assessment

The product has not been tested.

#### 12.6. Other adverse effects

No information available.

#### Further information

Avoid release to the environment.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Advice on disposal

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Wash with plenty of water. Completely emptied packages can be recycled.

### SECTION 14: Transport information

#### Land transport (ADR/RID)

##### 14.1. UN number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Inland waterways transport (ADN)

##### 14.1. UN number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Marine transport (IMDG)

##### 14.1. UN number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Micro Ortho Plus

Print date: 19.04.2017

Product code: 10-500\_EN

Page 6 of 6

#### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Air transport (ICAO)

##### 14.1. UN number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### 14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

#### 14.6. Special precautions for user

No information available.

#### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

not applicable

### SECTION 15: Regulatory information

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

##### EU regulatory information

##### Additional information

To follow: 850/2004/EC , 79/117/EEC , 689/2008/EC

##### National regulatory information

##### Employment restrictions:

Observe restrictions to employment for juvenils according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers.

##### Water contaminating class (D):

1 - slightly water contaminating

#### 15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

### SECTION 16: Other information

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of

Chemicals EINECS: European Inventory of Existing Commercial Chemical

Substances ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

#### Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)